

### ☒37. Synthetic Peptide - High Speed Screening

Sample: Ac-Lys-Tyr-Ala-Leu-Lys-Ala-Leu-Lys-Gly-Leu-Lys-acid  
Column: PLRP-S 100Å 3µm, 50 × 4.6mm (1512-1300)  
Eluent A: 0.1% TFA in 1% ACN  
Eluent B: 0.1% TFA in 50% ACN  
Gradient: 40 - 80% B in 5 mins  
Flow Rate: 1.0mL/min  
Detector: UV, 220nm

