

68. Typing Recombinant Monoclonal Antibody

Column: PL-SCX 1000Å 8µm, 150 × 4.6mm (1545-3802)
Eluent A: 10mM MES, pH 6.0
Eluent B: 0.2M NaCl, 10mM MES, pH 6.0
Gradient: 5% B for 5 min followed by
linear 5 - 85% B over 40 min
Flow Rate: 1.0ml/min
Detector: UV, 280nm

Reprinted from Journal of
Pharmaceutical & Biomedical
Analysis, Vol 16,
K G Moorhouse, W Nashabeh,
J Deveney, N S Bjork, M G
Mulkerrin and T Ryskamp,
*Validation of an HPLC method
for the analysis of the charge
heterogeneity of
the recombinant monoclonal
antibody IDEC-C2B8 after
papain digestion*, 1997,
with permission from Elsevier
Science.

