

☒ 17. Active Components from *Ginkgo biloba* Leaves

Column: PLRP-S 100Å 5µm, 250 × 4.6mm (1512-5500)

Eluent A: Water

Eluent B: ACN

Gradient: 5-25% B in 15 mins

Eluent C: THF

Gradient: 0-25% C in 25 mins

Flow Rate: 1.0ml/min

Detector: PL-ELS 2100

(neb=50°C, evap=70°C, gas=1.4 SLM)

KEY

1. Bilobalide

2. Ginkgolide A

3. Ginkgolide B

