

## ☒ 10. Sports Drink Analysis

Sample: High energy orange flavor,  
non-carbonated sports drink

Column: **PL Hi-Plex Na,**  
**300x7.7mm**

Eluent: Water

Temp: 80°C

Flow Rate: 0.3 ml/min

Detector: RI

POLYSACCHARIDES:

KEY

1. Dp8
2. Dp7
3. Dp6
4. Dp5
5. Dp4
6. Dp3
7. Dp2
8. Dp1(Glucose)
9. Fructose

